

Densfil® DensHeat® Oven III Obturator Oven

3-dimensional heat
for a 3-dimensional obturation



- 3-dimensional heat
- Indicator lights and audible alerts
- Slow release obturator holders
- Cleaning mode

Densfil® DensHeat® Oven III Obturator Oven

3-dimensional heat for a 3-dimensional obturation

3-dimensional heat

Obturator holders are heated in seconds in single and dual heating modes. 3-Dimensional heat is evenly maintained in the heating chamber until the obturator is ready to use, making it ideal for multi-canal procedures.

Indicator lights and audible alerts

The oven uses both LED indicator lights and audible alerts to indicate status changes.



Slow release obturator holders

Obturator holders are designed to allow stable vertical entry and exit of the heating chamber.

Cleaning mode

The cleaning mode raises the temperature of the oven and softens gutta-percha inside the heating chamber for removal.

- The oven heats up obturators in 20 to 49 seconds, depending on the size of the obturator
- A built-in timer audibly signals when the filling material reaches the right temperature
- The obturator is kept at the right temperature for 90 seconds.

Item Number	Contents
DENSHEATIII	Densfil® DensHeat® Oven III Cleaning brush
DENSHEATIIISYS	1 Densfil® DensHeat® Oven III 20 Plastic Obturators (20-40 Asst.) 1 Size Verification Kit Refill (20-40 Asst.) 1 AH Plus Jet® Root Canal Sealer Starter Kit 1 Practice Kit 1 Densfil® Manual