

Increased Taper. Incredible Value.



.04 & .06 Tapered Gutta Percha and Absorbent Points.

Designed to match today's tapered shaping instruments, Maillefer's .04 and .06 Gutta Percha and Absorbent Points deliver the precision you need for efficient, superior results. And because they are packaged to give you more for your money, you get an exceptional deal.

Absorbent Points

- Designed to fit .04 or .06 prepared canal
- Sold in packages of 180 points
- Hand-rolled
- Available in ISO sizes 15-40

Gutta Percha Points

- Allows for a single cone obturation
- Fits securely into .04 or .06 instrumented canal
- Sold in packages of 120 points
- Hand-rolled
- Available in color-coded ISO sizes 15-40

Order Today.

800-924-7393 • 800-924-7389 (fax) • www.maillefer.dentsply.com

DENTSPLY
MAILLEFER

In Endodontics, Quality Has A Name.

MAILLEFER

Maillefer .04/.06 Tapered

Gutta Percha

and Absorbent Points

Gutta Percha .04 Taper

670948	#15	White
670949	#20	Yellow
670950	#25	Red
670951	#30	Blue
670952	#35	Green
670953	#40	Black
670954	#15-40	Assorted

.04

Gutta Percha .06 Taper

670955	#15	White
670956	#20	Yellow
670957	#25	Red
670958	#30	Blue
670959	#35	Green
670960	#40	Black
670961	#15-40	Assorted

.06

Absorbent Paper Points .04 Taper

670186	#15	White
670187	#20	Yellow
670188	#25	Red
670189	#30	Blue
670190	#35	Green
670191	#40	Black
670192	#15-40	Assorted

.04

Absorbent Paper Points .06 Taper

670193	#15	White
670194	#20	Yellow
670195	#25	Red
670196	#30	Blue
670197	#35	Green
670198	#40	Black
670199	#15-40	Assorted

.06

Order Today.

800-924-7393 • 800-924-7389 (fax) • www.maillefer.dentsply.com